

Environmental Sample Processor Mission Scripting



5/22/14 Brent Roman brent@mbari.org

Mission Scripts and Phases

- Top Level Commands for a deployment

 Often omitted for lab work
- Usually contains a mission method
 - Specifies the starting tube number
 - Optionally specifies Mission End Time
 - Contains any number of mission phases
 - Each having a start time
 - with optional trigger conditons
 - One or more protocols run per phase
 - The ESP sleeps between phases
 - Contextual sensors continue being polled



Protocols

- Protocol scripts do the real work of microbiological assays
 - Many canned scripts available:
 - HAB = Harmful Agal Bloom
 - BAC = Bacterial Assay
 - LARV = Larval Assay
 - WCR = Whole Cell Archival
 - DA = Domoic Acid Assay
 - HABDA = combined HAB and DA assay
 - STX = Saxitoxin Assay
 - All have parameters you may modify to suit your needs
 - With default values so you needn't specify everything
 - You may also create new protocols using the existing protocols as a guide:
 - STX was created just last year as a variant of DA



Example "3peat" QC Mission

mission startTube: 2, endTime: "6AM 12/18/12" do

```
at "12:40:00 12/14/12" do
habda {noKill}
end
```

```
at "12:40:00 12/15/12" do
habda {noKill}
end
```

at "04:00:00 12/17/12" do habda end



Copyright MBARI 2014

It's Ruby all the way down

- Commands, Missions, Scripts, Protocols, Configuration Files
 - All are written in version 1.8 of the Ruby scripting language
- Learn a little Ruby
 - * Rote memorization fails when something goes \underline{W} rong
 - Standard on Mac OS, easily installed everywhere else.
- A gentle tutorial:
 - https://pine.fm/LearnToProgram/
- The bible:
 - http://pragprog.com/book/ruby/programming-ruby
- More (free) choices to suit your learning style:
 - http://ruby.about.com/od/tutorialsontheweb/tp/10waysfree.htm

