



Physical State

5/24/22 Brent Roman brent@mbari.org



ESP Physical State

- Canister environmental sensors
- Power Switches
- All moving actuators:
 - Rotary Valves
 - Solenoid Valves
 - Syringes
 - Clamps
 - Carousel
 - Elevator
 - Elbow
 - Gripper
- Puck Heaters



Gateway Canister State

- Quietly logged once every

- 10min for 2G, 2min for 3G

-> `Gate.queryCan` #or simply -> `can`

```
Can@22:40:11, 24.2C, 67% humidity, 14.2psia, 13.684V, 0.256A,  
187.834Ah, 3.50W
```

- Ah is only available with newest firmware

-> `Gate.queryCan` takes an immediate reading

- `can` is an alias for `Gate.queryCan`

-> `Gate.can` uses the most recent reading

-> `Gate.canPollInterval=30.seconds`

-> `Gate.canPollInterval=0` #no polling



Canister State details

-> `Can.unit` #all available can sensors

```
{avgCurrent: {format: "%.3fA avg"},  
  batteryUsed: {format: "%.3fAh"},  
  current: {format: "%.3fA"},  
  humidity: {format: "%d%% humidity"},  
  pressure: {format: "%.1fpsia"},  
  temperature: {format: "%.1fC"},  
  voltage: {format: "%.3fV"},  
  waterAlarm: {format: "%d%% Wet!"},  
  threshold: 1.5}}
```

-> `Can.missionDuration` #duration power has been on

```
5 weeks, 3 days, 10:19:32.40625
```

-> `c=can; (c.voltage*c.batteryUsed).round` #~Wh used



Gateway Power Switches

-> `Power` #state of all power switches

```
{analytic1: false,  
  analytic2: false,  
  analytic3: false,  
  camera: false,  
  core: false,  
  raw: false,  
  sampler: false}
```

-> `Power.sampler :on` #turn sampler on

-> `Power.on :sampler` #turn sampler on

-> `Power.off :sampler, :camera` #turn off sampler and camera

-> `Power.offSince[:sampler]` #when sampler was turned off or nil

-> `Power.onSince[:core]` #when core was powered on or nil

-> `Power.sampler` #True if sampler is on



Linear Actuator Class Hierarchy

- Slide => Lowest level at which end-users interface with hardware
 - Think of a slide trombone with named positions for arbitrary “notes”
 - Forearm, Elbow, Carousel
- Clamp => inherits from Slide
 - Adds closed?, open?, and closeAndVerifyPuckPresence
- Scale => inherits from Slide
 - Adds a linear, numeric scale (in standard units) to Slides
 - Elevator #unit = puck height
- Syringe => inherits from Scale
 - adds pull, push, fill, empty volume methods
 - Collection, Processing, Sampler, Analytical syringes
- Thermal => inherits from Scale
 - replaces motor with a heater
- Errors come from microcontrollers, which are all managed by Slide class
 - This is why Linear Actuators report most errors as “Slide::Error”



Other Actuator Classes

- Shaft is a rotary actuator
 - all of which just happen to spin a rotary valves
- Gripper is a two state actuator
 - with binary position sensing
 - May control a robotic hand or a motorized valve
- Solenoid is a two or three state solenoid actuator
 - without position feedback sensing
- Valve is a Solenoid used to control fluid flow
- Valve::Manifold is an array of Valve
 - sharing a common fluid path



Using AxisMaps

- An [AxisMap](#) maps all raw counts to corresponding position names
- They are typically accessed via their associated Axis or Positions:
 - `axis.legend` => the AxisMap as a Hash
 - `axis.list` => list of all names without raw positions
 - `axis.labels` => list of only the position labels – omitting aliases
 - `axis.maxPosition` => position mapped to greatest raw counts
 - `axis.minPosition` => position mapped to least raw counts
 - `axis.advance` => move to position with next higher raw counts
 - `axis.retard` => move to position with next lower raw counts
 - `axis.at?(position)` => true if axis is at (or near) specified position
 - `axis.near?(position)` => true if axis is at or near position
 - `axis.between?(pos1,pos2)` => true if axis is (nearly) between
 - `axis.rawId(rawCount)` => position nearest rawCount (reverse map)
 - `position.advance(detents)` => position with next higher raw counts
 - `position.retard(detents)` => position with next lower raw counts
 - `position.near?(position)` => true if positions very near each other



Axis Map Example

- Hardware counts \Leftrightarrow names and aliases

-> `Forearm.legend`

```
{-12793 => 3, -12782 => 5, -12726 => 7, -12711 => 4,  
-12706 => 2, -12699 => 1, -12678 => 6,  
-12501 => [CC, :collection, Collection],  
-12418 => [PC, :processing, Processing],  
-12342 => [:garage, FlushPuck, FlushPuck::Garage],  
0 => :home,  
2800 => ["retracted", :retract, :clear]}
```

- Defined in `configure.rb` as:

```
Forearm.detents 0=>:home, 2800=>"retracted",  
-12501=>CC, -12418=>PC, -12342=>:garage,  
-12699=>1, -12706=>2, -12793=>3, -12711=>4, -12782=>5,  
-12678=>6, -12726=>7
```



Solenoid Valves

- Solenoid valve state is either open, closed, or unknown
 - Intake and Exhaust external solenoid valves
 - Unipoler, normally closed
 - Draw a lot of power while opened
 - Solenoid valves in Manifolds
 - Bipoler, latching
 - Latching Valve state is initially unknown
 - Draw power in short pulses when changing state
- `Intake.to :open` #holds Intake open
- `Intake.open` #ditto
- `Solenoid` #shows Intake and Exhaust state



Solenoid Valves Plumbing

- Each ESP 2G Dwarf microcontroller
 - drives 8 solenoids numbered 0..7
- To show how a solenoid is connected:
 - `Intake.wiring` #how is Intake wired?
 - `sampler[6]` #it's 2nd to last on sampler dwarf



Valve::Manifold

- Each composed of a series of Solenoid::Valve plus an endName
 - Valve::Manifold state = either the name
 - of its first opened Solenoid::Valve
 - or its endName, if no Solenoid::Valve is open
 - Individual valves in the manifold can be accessed
 - `CSR.series[0].open` #open 1st CSR valve
 - `CSR.wiring` #shows wiring connections
- ```
[:lysis <=> collection[0], #these need not be in order
 :diluent <=> collection[1],
 :RNAlater <=> collection[2],
 :mfbkill <=> collection[3],
 :kill <=> collection[4],
 :flush <=> collection[5]]
```



# Valve::Manifold Configuration

- Collection Series configuration is often machine specific
- Example for ESPwaldo, defined in its configure.rb file:

```
:CSR.denotes Valve::Manifold :Collection, [
 Valve.reagent(:lysis, CollectionValves,0),
 Valve.reagent(:diluent, CollectionValves,1),
 Valve.reagent([:RNAlater, :rnal], CollectionValves,2),
 Valve.reagent(:mfbkill, CollectionValves,3),
 Valve.reagent(:kill, CollectionValves,4),
 Valve.reagent(:flush, CollectionValves,5)], :air
```

- :air is the manifold's "endName" representing
  - its state when all its series valves are closed
- :rnal is an alias for :RNAlater



# aliases and labels

- Each named position has exactly one label
  - in addition, it may have any number of aliases
- A position is output by its label
  - but may be input by its label or any corresponding alias

```
-> CSR.aliases
```

```
{rnal: :RNAlater}
```

```
-> CSR.alias :clean=>:flush #clean now an alias for flush
```

```
-> CSR.relabel :atmosphere, :air
```

```
-> CSR.to :air
```

```
Collection Valve::Manifold selects atmosphere
```

- rather than air!



# Rotary Valves

- Each rotary valve is controlled by its Shaft position
  - Shaft is the class of all rotary valves
- Raw shaft state is an angle from 0..511
  - A number of these may be defined as named positions
- Rotation direction to the new goal position
  - may move over raw position 0
  - may be specified
    - to avoid moving over another position
- A goal position may be specified
  - as a named position
  - exactly between 2 named positions
  - as a raw offset from one of the above



# Shaft Configuration

- Each 2G Dwarf may control up to 4 rotary valves numbered 0..3

-> `PTV.wiring` #displays wiring information

`processing[2]` #wired to 3rd channel of processing dwarf

-> `PTV.legend` #displays position map

`{64 => [PRV, 1, #positon labeled PRV [with alias 1] is raw angle 64`

`192 => [:PRVmixing, 2],`

`320 => [:mixing, 3],`

`448 => [:puck, 4] }`





# Rotating Shafts

- > `PTV.to :mixing #rotates to mixing position fastest way 'round`
- > `PTV.select :mixing, avoiding: PRV`
  - #rotates to mixing position
  - in a direction that avoids rotating by the PRV position
- > `PTV.select :mixing, via: PRV #rotates opposite way!`
  - #rotates to mixing position via the PRV position
- > `PTV.dialBetween :mixing, :puck`
  - #rotates to between mixing & puck position fastest way 'round
- > `PTV.dialBetween :mixing, :puck, avoiding: PRV`
  - #rotates to between mixing & puck position, avoiding PRV
- > `PTV.rawAngle #return the raw angular position of the shaft`  
320
- > `PTV.at? :mixing #true if PTV is at, or very near, mixing`  
true



# Gripper Characteristics

- Two position actuator with minimal position feedback
- In one of two states
  - or transitioning between those states
- The two states are named when configured
  - they may not have aliases or labels
- Examples are:
  - The ESP 2G Hand was original Gripper actuator
  - Some ESP 2G external Sample valves are controlled as Grippers
  - 3G External rotary valves are all Grippers

# Gripper Use

- > `Hand.close` #closes the Hand Gripper
- > `Hand.open` #opens it
- > `Hand.state` #:open, :closed, or :unknown
- > `Hand.open?` #true if Hand is open

true

- Only on NOAA GLERL 2G ESPs...

- > `Sample.deep` #moves Sample valve to its :deep position
- > `Sample.shallow` #moves Sample valve to its :shallow position
- -> `Sample.state` #:shallow, :deep, or :unknown
- -> `Sample.deep?` #true if Sample valve is in its :deep state

false



# Slide Characteristics

- Linear Actuator having precise position feedback
- Raw positions are in hardware counts
  - Each Slide's count units may be different
  - Most Slides require physical 'homing'
    - to calibrate their counts position sensor
    - A Slide that is not yet homed is 'lost'
  - are mapped to names with associated AxisMap
- Each Slide includes a set of named configurations
  - that set motor limits and velocity profile
  - only one such configuration is active at a time

# Basic Slide Operations

- The Slide is the “base class” for linear actuator axes
- `slide.configure cfg =>` forces configuration object `cfg` to dwarf
- `slide.reconfigure cfg =>` sends `cfg` only if changed from last
- `slide.in(cfg) {block} =>` execute block in configuration `cfg`
- `slide.position =>` return the slide's current position
- `slide.goal =>` return the slide's current goal position
- `slide.jog counts =>` move specified # of raw encoder counts
- `slide.seek goal =>` move to specified goal position
  - Without updating servo's configuration
- `slide.to goal, config =>` move to specified goal position
  - Updating servo's configuration if appropriate
- `slide.hold =>` hold the current position
- `slide.coast =>` turn off the servo
- `slide.force =>` apply constant “force” (`slide.force 0 = slide .coast`)
- `slide.stop =>` brake to a stop as fast as possible
- `slide.log(decimator) {block} =>` log slide status while doing block
- `slide.status =>` return current slide servo status object



# How Scales Differ from Slides

- Scales inherit all the operations of Slide, adding:
  - Linear mapping of logical “amounts” or “units” to raw counts
    - $\text{rawCount} = \text{scale.countsPerUnit} * \text{amount} + \text{zero}$
    - zero is simply the rawCount value at 0 amount
    - `scale.zero` => -12580 #example case
    - `scale.gain` => `scale.countsPerUnit` => 32498.0
- AxisMap associated with a Scale:
  - Must contain at least two positions whose labels are numeric
    - [there should be only two numerically labeled positions]
    - These positions project the scale's linear mapping onto counts



# Scale::Skew objects

- A Scale::Skew is a generic, linear mapping object
  - represents
$$y = mx + b$$
  - `scale.skew` => `-12580.000+32498*counts`
  - `scale.skew.gain` => `32498.0`, `scale.skew.bias` => `-12580`
  - `scale.skew.apply(2)` => `52416.0` # == `2*32498 - 12580`
    - solves for y (engineering units) given x (counts)
  - `scale.skew.reverse(scale.skew.apply(x))` => `x`
    - solves for x (counts) given y (engineering units)
  - `Skew.bestFit(counts, units)` => skew that best fits data
  - `Skew.interpolate()` => interpolates among array of skews
- Scale::Skews are also used to calibrate Thermal pads!



# How Syringes differ from Scales

- A syringe is merely a scale with volumetric units
- volume is defined as an alias for amount
- Similarly for maxVolume and minVolume
- fill method moves to the syringe's maxPosition
- empty method moves to the syringe's minPosition



# 2G ESP Dwarf DC Motor Servos

- Two identical servo channels
- 64hz sampling timebase (sample rate typically 32hz)
- Each Channel's Inputs:
  - Quadrature incremental encoder
    - (A and B 90 degrees out of phase)
  - Home flag (typically a hall effect sensor)
  - Optional threshold sensor
  - Forward and Reverse limit switches
  - One General Purpose digital input bit (for gripper)
- Each Channel Outputs:
  - PWM -100% to 100% (15 kHz with 1% resolution)
  - One General Purpose digital output bit



# Configuration Object Details

- :samplePeriod = number of 64hz timebase tics per sample tic
  - Default value = 2 (Typically 1 or 2)
- :encoder, :threshold, :home sensor power / polarity
  - Default value = :off (may be :positive or :negative)
- :homeDirection = :forward or :reverse
  - Default value = :reverse
  - :reverse moves negative if home flag inactive
- :brake = short motor terminals on servo error (:false or :true)
  - Default value = true
- :debug = output servo state at sample rate while seeking goal
  - Default value = false



# Control Gains and Factors

- PID :gain struct with members P, I, and D
  - Default values for each are 0
  - Servo will not operate until at least one is non-zero
  - Effective value of P and D is divided by 4096
  - I is effectively divided by 16384
- :friction compensation gain
  - $\text{cmdVel} * \text{friction} / 4096$  added to PWM output
  - $\text{cmdVel} = \text{Commanded velocity}$
- :stiction compensation factor
  - If negative  $\text{cmdVel}$ , subtract  $\text{stiction}/2$  from PWM
  - If positive  $\text{cmdVel}$ , add  $\text{stiction}/2$  to PWM

# Trajectory Generator (1 of 2)

- :acceleration & :deceleration in counts/tic/tic
  - Default values for each are 0, normally positive
  - Specify negative acceleration to disable “softstart”
  - Zero :deceleration implies  $\text{deceleration} = \text{abs}(\text{acceleration})$
- :maxSpeed = plateau velocity in counts/tic
  - Temporarily reduced when PWM limits reached to prevent trajectory errors due to low battery voltage
- :minSpeed = slowest acceptable progress rate (counts/sec)
  - Speed error if maxSpeed reduced below minSpeed
- :maxSettling = max tics to allow to servo to settle at goal
  - Default 0, typically 2 – 3 seconds worth of tics
  - Just ensures that position error not returned too early



# Trajectory Generator (2 of 2)

- `:stopWindow` determines how nearly goal should be reached
  - Specified in encoder counts (16 bit limit max)
  - Temporarily increased each time goal is passed
  - Special Value false indicates no (more) reseek allowed
  - Defaults to Special Value `:deceleration` = deceleration rate
  - Also accepts value `:acceleration`
- `:hunt` determines whether to adjust setpoint after goal reached
  - Defaults to false, set true to “fight” to hold exact position at goal
  - Setpoint is *never* adjusted if position within `stopWindow`
- `:thresholdOffset` determines how far from threshold to stop when reached
  - Defaults to 0 encoder counts
  - When threshold reached before goal,  $goal = position + thresholdOffset$
  - Used to position top of puck stack with respect to ESP's top plate



# Core Limits

- :maxPWM & :minPWM
  - Max must be  $\geq$  min, but each may be negative or positive
  - Constrains servo output, but does not constrain “force” command
  - Effective maxSpeed is reduced when servo reaches these PWM limits
- :maxPositionErr determines absolute maximum tolerable servo error in different contexts:
  - SeekErr if stopWindow grows too large due to repeatedly missing goal
  - TrajectoryErr if position becomes too far from setpoint while transiting
  - PositionErr if position moves too far from goal after arrival
- :maxCurrent determines maximum allowable motor current
  - In milliamps
  - Should never be set  $> 2000\text{mA}$



# Pressure Limits

- :maxInPress, :maxOutPress, :minInPress, :minOutPress
  - 0 to 4095 ADC counts
  - Maximum/Minimum tolerated Intake and Outlet pressures
  - Constraint disabled if corresponding max == min
  - All default to 0
- :maxDeltaPress & :minDeltaPress -- (-4095 to 4095)ADC counts
  - Maximum/Minimum tolerated pressure difference
  - Constraint disabled if set to special value: false
  - All default to false (there is no corresponding value true)
- Generic “Pressure Error” results if any of the above are violated
  - One must check status to determine the exact problem



# Pressure Servo Configuration

- `:inputDeltaPress` determines if pressure delta is sensed or derived
  - True to input the difference from ADC 7
  - False to derive it as (intake – outlet) pressure
  - Defaults to false
- `:pressBias` is subtracted from delta pressure before use
  - In servo or limit check
  - Defaults to 0
- `:pressGain` is the proportional gain of the pressure servo
  - Scaled like P and D, `pressGain` is \*4096
  - Reduces acceleration from that normally determined by the trajectory generator.
  - Never causes command velocity to fall below `minSpeed`





# Default Processing Syringe I2C::Servo::Configuration

```
:PS.denotes Syringe "Processing Syringe",
 :processing, 0,
 :encoder=>:negative, :home=>:negative,
homeDirection:false,
maxPositionErr:65,
gain:PIDgain(3500, 3000, 1300),
friction:170,
maxSpeed:100, minSpeed:30,
acceleration:5,
maxCurrent:120,
maxSettling:3*32
```

`PS.maxDuration=160` #how long can a move take?

- *excerpted from shallow/preconfig.rb*



# Alternative Processing Syringe Servo Configurations

- Based on default configuration shown on previous slide:  
PS.defCfg :fast, maxSpeed:300  
PS.defCfg :slow, maxSpeed:50  
PS.defCfg :slow1, maxSpeed:10, minSpeed:2, acceleration:2
- Based on :slow1 configuration (defined above):  
PS.defCfg :slow2, :like=>:slow1, minSpeed:5
- *excerpted from shallow/postconfig.rb*



# Switching Servo Configurations

### The *easy* (and correct way) ###

```
PS.to PS.maxVolume/2, :slow1 #half full (or is it empty?)
```

- Only changes the configuration if necessary
- Don't use `.seek` unless sure the config already loaded on dwarf.

### The hard (and also correct way) ###

```
PS.in :slow1 do
```

```
 PS.to PS.maxVolume/2
```

```
 PS.empty #this is still in PSslow1
```

```
end
```

```
PS.fill #old configuration restored (likely PSconfig)
```

- `slide.in {block}` constructs may be nested arbitrarily deep



# I2C::Servo::Status Objects

- :enabled = true if servo control is active
- :pastRLS, :pastFLS, :pastThreshold, :home
  - True if corresponding switch is closed
- :position = 32-bit signed offset from home position
- :velocity = 16-bit signed in encoder counts/tics
- :current = signed milliamps
  - Always agrees with sign of PWM status below
- :PWM = signed percent PWM duty cycle
- :err = 16-bit signed (setpoint – position)
- :voltage = raw motor voltage (in volts)
  - This is the *only* floating point value



# Servo Pressure Status

- Recall that pressure may be a proxy for any arbitrary voltage input
- :inPress = intake pressure in raw ADC counts (0-4095)
- :outPress = outlet pressure in ADC counts
- :deltaPress = delta pressure in ADC counts
  - This is always ADC channel 7
  - It is *not* affected by the :inputDeltaPress configuration flag

# Plotting Slide Servo Trajectories

- Add ssh key to workstation's `authorized_keys` file
  - So that the ESP host can run commands without password prompts
    - Test from Linux shell prompt, on ESP host, by invoking:  
`$ ssh workstation ls`
    - *This is a security breach. Remove key when done if it worries you.*
- Edit `remotePlot` method `utils/plot.rb` as necessary
  - To change the workstation name (and possibly the display number)  
-> `require 'plot' #only once per session`
- To produce each new plot window:  
-> `remotePlot slide.log {blockOfCodeExercisingSlide}`
- e.g. plotting default status fields of position, velocity and current:  
-> `plot CC.log {CC.to :closed}`
- e.g. plotting `:current`, `:voltage`, `:pwm`, and `:err`  
-> `plot(CC.log {CC.to :closed}, :current, :voltage, :pwm, :err)`

