



ESP Fluid Tracking

5/29/22 Brent Roman brent@mbari.org



What is Fluid Tracking?

- The 2G ESP can follow fluid movements
 - between Syringes and containers
 - into and out of clamped pucks
 - into and out of external Intake and Exhaust
- Fluid tracking works in real and simulated modes
 - can be configured to stop operation
 - just before a reagent bag is “pulled flat”
- Fluid tracking slows simulation by about 20%

Using Fluid Tracking

- > `require 'fluid'` #enables fluid tracking
 - > `Tank.clear` #empties all virtual containers
 - Run ESP normally... later...
 - > `Tank` #displays long table summarizing fluid use
- ```
#Fluid::Vessel# Net Volume (Lowest Volume)
```
- `Fluid::Vessel` = name of the fluid container
  - `Net Volume` = change in volume of fluid since `Tank.clear`
  - `Lowest Volume` = least volume remaining in container
    - only tracked for containers where fluid is pumped both ways
  - first column of each line is a special code
    - '=' denotes closed fluid container (eg. :flush)
    - '~' denotes container that is vented to atmosphere (eg. ~waste)
    - ':' denotes fluid added to or removed from above container

# Limiting Reagents' Supply

```
-> require 'fluid' #enables fluid tracking
```

```
-> Tank.clear #empty all virtual containers
```

```
-> Tank.supply flushCSR:12, flushPSR:5
```

```
 #supply only 12ml flushCSR, 50ml flushPSR
```

```
-> bac 100 #start the usual bac protocol
```

```
...
```

```
Pulling 5ml of CSR.flush
```

```
CSR.select :flush
```

```
CS.seek 5ml
```

```
Fluid::Depleted in simfast -- flushCSR pulled dry by Collection Syringe
```

```
-> Tank.supply flushCSR:20 #hang a new virtual flush bag
```

```
-> Thread.resume #resume from Fluid::Depleted error
```



# Banking Virtual Reagents

```
-> require 'fluid' #enables fluid tracking
-> Tank.clear #empty all virtual containers
-> Tank.bank flushCSR:12, flushPSR:5
 #deposit balance of 12ml flushCSR, 50ml flushPSR
-> bac 100 #start the usual bac protocol
```

...

- Banked reagents are not depleted when they run out
  - instead, their balance simply goes negative!
- If you use Tank.supply and Tank.bank for different reagents
  - only those you last “supplied” can raise Fluid::Depleted errs
- Tank.clear zeros the balance of all reagents



# Logging Fluid Movement

- > `require 'fluid'` #enables fluid tracking
- > `Tank.clear` #empty all virtual containers
- > `Fluid.logTransfer true` #log fluid movements
- > `bac 100` #start the usual bac protocol

... Log excerpted during Puck Evac ...

```
<simfast> CBV.dial SSV
```

```
@23:29:16.16 <connect> CC-->SS 3.341ml @ 3.01bar
```

```
CS-->SS 0.792ml @ 3.58bar
```

```
<simfast> CTV.dial :air,avoiding: :intake
```

```
@23:29:17.66 CC-->air 10.072ml
```

```
SS-->CC 2.817ml @ 2.29bar
```

```
SS-->air 6.454ml
```

...

