

Environmental Sample Processor Mission Scripting



5/25/22 Brent Roman brent@mbari.org

Mission Scripts and Phases

- Top Level Commands for a deployment

 Often omitted for lab work
- Usually contains a mission method
 - Specifies the starting tube number
 - Optionally specifies Mission End Time
 - Contains any number of mission phases
 - Each having a start time
 - with optional trigger conditons
 - One or more protocols run per phase
 - The ESP sleeps between phases
 - Contextual sensors continue being polled



Protocols

- Protocol scripts do the real work of microbiological assays
 - Many canned scripts available:
 - hab = Harmful Agal Bloom
 - bac = Bacterial Assay
 - larv = Larval Assay
 - wcr = Whole Cell Archival
 - da = Domoic Acid Assay
 - habda = combined HAB and DA assay
 - stx = Saxitoxin Assay
 - All have parameters you may modify to suit your needs
 - With default values so you needn't specify everything
 - You may also create new protocols using the existing protocols as a guide



Example "3peat" QC Mission

mission startTube: 2, endTime: "6AM 12/18/12" do

```
at "12:40:00 12/14/12" do
habda {noKill}
end
```

```
at "12:40:00 12/15/12" do
habda {noKill}
end
```

```
at "04:00:00 12/17/12" do
habda
end
```



end

Example Mission (1 of 3)

Top level ESP Mission script for July 2017 Niagara deployment - niagara17jul11_V2.rb
\$daVol = 1000 #5.0um durapore

\$da[:wcrVol]=1000 #default to 1L wcr for da
\$startTime="+2 days" #relative to start of last phase – not its end time!

\$numPhases = 0 #Deployment count of phases run

def nextTime startTime=\$startTime

#if a relative time is specified, reference the start of the

#previous phase instead of the current time

Delay.wakeTime startTime, \$started

end

Copyright MBARI

```
nextStartTime = method :nextTime
```



Example Mission (2 of 3)

def withoutWCR [\$daVol, nil] end

def twoDAs logPhaseStart :twoDAs Sample.shallow; initialPurge; da withoutWCR Sample.deep; initialPurge; da withoutWCR end

fill! if ESP.simulation? #simulate full puck load



Example Mission (3 of 3)

mission :startTube=>2 do

```
at "4PM 7/17/2017" do
```

Sample.shallow; initialPurge; da withoutWCR

```
end
```

```
at "10AM 7/19" do
```

```
Sample.shallow; initialPurge; da withoutWCR
```

```
end
```

```
3.times {
```

at nextTime do

```
Sample.shallow; initialPurge; da withoutWCR end
```

```
}
```

```
3.times {
```

```
at nextTime '10AM' do
```

twoDAs

end

}

15.times {

at nextTime do

twoDAs

end

}

end Copyright MBARI 2022



Running Late

- When the phase start time has already past
 Delay::MaxLate determines whether
 - to throw a TooLate Exception, or
 - to continue with a running late warning
 - defaults to 1.hours

-> Delay.adjust :MaxLate, 2.hours

• To limit the longest delay to start next phase

-> Delay.adjust :MaxWait, 2.weeks

- MaxWait defaults to nil
 - which allows mission phases to wait forever



Alternate Framework for daily repeating phases

- Schedule is an array of daily start times

 rather than being coded in Ruby mission script
- Assumes next day is like the present one
- Assumes startTube already initialized
- Does not start and stop contextual sensors
- Every phase must process the same number of pucks
- Easy to skip weeks or months
 - combersome if every day is different
- It's easy alter schedule from an espclient
 - without aborting the running mission
- Mission runs until all but a set number of pucks processed
 - conserves these for controls



Example Daily Mission

```
require 'daily' #use alternate "daily" mission framework
Daily.mission 8, 4, [ #reserve 8 pucks, process 4 at a time
    "2/7/23".daily("9AM", "3PM", "11PM"),
    "2/8".daily("rise + 3", "set + 2:15:00"),
    "2/12".daily,
    "2/15".daily("9:14", "19:59")
] do |startTime, remainingPucks|
    at startTime do
        initialPurge; bac
    end
```

end

- In English: while more than 8 pucks remain,
 - 3 bac phases each day on 2/7/23
 - 2 bac phases each day from 2/8 through 2/11
 - sleep from 2/12 through 2/14 #vacation :-)
 - 2 bac phases each day from 2/15

