



2G ESP Water Sampling

5/25/22 Brent Roman brent@mbari.org



Syringe Pump Driven Sampling

- The 2G ESP uses its 25ml Sampling Syringe
 - to pull water through filter in Collection stage
 - while monitoring syringe pressure
- Syringe stops when pressure in it falls too low
 - until pressure equalizes
 - resumes at 90% of the previous pumping rate
- Different Samplers optimized for each filter type
 - [BACsampler](#), [HABsampler](#), [DAsampler](#), etc.
 - see `utils/sampler.rb`, `utils/gensampler.rb`,
 - and `protocol/2G/[BAC,HAB,DA].rb`

2G Sampler Parameters (1 of 5)

Defined in utils/gensampler.rb:

- **vacuum** defaults to 10..13 [PSI]
 - target vacuum range while pulling through filter
- **minRate** defaults to 0.2 [ml/second]
 - Filter is clogged when filtering at minRate
 - generates vacuum > vacuum.max
- **avgRate** defaults to 0.2 [ml/second]
 - Constant filtering rate (used only during simulation)

2G Sampler Parameters (2 of 5)

Defined in utils/sampler.rb:

- **bubblePt** defaults to nil [PSI]
 - Filter's characteristic bubble point (*must be spec'd*)
- **evacVolume** defaults to 10 [ml]
 - How far to pull back SS to evacuate puck
- **endVacuum** defaults to 3 [PSI]
 - Vacuum must fall below this before starting next stroke
- **maxDelay** defaults to 40.seconds [seconds]
 - How long to wait for endVacuum to be reached
 - Filter Clogged if vacuum still > endVacuum after maxDelay

2G Sampler Parameters (3 of 5)

Defined in utils/sampler.rb:

- **backoff** defaults to 0.9 [seconds/seconds]
 - new filtering rate after restarted due to high vacuum
- **evacDelay** defaults to 40.seconds [seconds]
 - how long maintain high vacuum during puck evac
- **numberOfEvacs** defaults to 1
 - number of times to repeat puck evacuation
- **maxPumpVolume** defaults to 40 [ml]
 - max vol of air pumped to reach bubblePt during evac

2G Sampler Parameters (4 of 5)

Defined in utils/sampler.rb:

- **equalized** defaults to 1.0 [PSI]
 - Intake blocked if pressure will not equalize to within this psi
- **minPSI** defaults to 2.5 [PSIa]
 - Minimum sampling pressure (to avoid boiling sample)
- **maxResidue** defaults to 0.1 [ml]
 - don't take another stroke if remaining vol < maxResidue
- **primeVolume** defaults to 7 [ml]
 - volume of initial prime stroke
- **maxPumpVolume** defaults to 40 [ml]
 - max vol of air pumped to reach bubblePt during evac



2G Sampler Parameters (5 of 5)

Defined in utils/sampler.rb:

- **meterInterval** defaults to 30.seconds [seconds]
 - Flow meter reporting interval
- **meterDelta** defaults to .1 [Liter/min]
 - report flow changes only if > meterDelta
- **exhaleConfig** defaults to nil [I2C::Servo::Configuration]
 - SS configuration for exhalation

Sampling in the Lab

```
-> mySampler = Sampler.new bubblePt:26
... verify parameters are to your liking ...
-> mySampler.maxDelay=60.seconds #increase maxDelay
-> CC.loadPuck :sh1 #load puck from startTube
-> mySampler.prime #take the optional prime stroke
-> mySampler.sample 55 #sample up to 55ml more
-> mySampler.avgRate=2 #simulate faster!
-> mySampler.sample 125 #sample up to 125ml more
-> mySampler.puckEvac #dry the puck
-> CC.unloadPuck #put spent puck into discard tube
```


collectSample

```
-> CC.loadPuck(:sh1) {collectSample(mySampler).puckEvac}
```

- Loads sh1 puck from startTube
- primes and samples default goal volume with mySampler
- dries the puck
- discards it into startTube-1
- replace mySampler with any of the predefined samplers
 - **BACsampler**, **HABsampler**, **DAsampler**, etc.

Sampler Pressure Sensing

- > `Sampler.intakePSI` #abs pressure above the filter
- > `Sampler.outletPSI` #abs pressure below the filter
- > `Sampler.PSI` #[abs pressure above, abs pressure below]
- > `Sampler.calibrate` #tare pressure sensors
- A dry puck must be clamped in the collection stage
- Must tare after a new puck is clamped

- For hardware debugging and initial calibration:
 - > `CS.status.outPress` #intake pressure in raw ADC counts
 - > `SS.status.outPress` #outlet pressure in raw ADC counts
- Raw pressure ADC counts range 0..4095

External Sampling Aids

- External (Dumb) Pumps
 - Switched Lift Pumps
- Intake selector valves
 - Choose between shallow and deep sample intakes
- Passive Flow Meters
 - Verify water flowing through sample intake loop
- Active Flow Controllers
 - Combine meter with pump controller to regulate flow
 - Also controls an external intake selector valve

Configuring External Pumps

- Switch via ESP's :sampler external power switch
- Sampler::PumpPrimeDelay
 - defines how long pump runs to prime external loop
 - before ESP intake opens
 - typically defined in machine's configure.rb as:

```
class Sampler
  const :PumpPrimeDelay, 3.minutes
end
```

Intake Selector Valve

- External Valve selects between deep and shallow intakes
 - Only installed on GLERL ESPs
 - Connected as a Gripper
 - to the Sampler Dwarf microcontroller
 - Example config at end of ESPnessie's configure.rb:

```
require 'depthsel'
```

```
:Sample.denotes DepthSelectorValve "Depth Selector Valve", :sampler, 1
```

```
Sample.maxDuration=30 #allows up to 30 seconds for valve to move
```

- Valve must be explicitly positioned before sampling begun
 - > `Sample.shallow` #select shallow intake
 - > `Sample.deep` #select deep intake



Passive Flow Meters

- Two types in use:
 - Analog input (by WHOI)
 - may function as a simple flow presence/absence switch
 - or uncalibrated % of full flow
 - in ESP's Can data
 - Serial input (by NOAA NWFSC/APL)
 - on ESP's external sampler RS-232 serial port
 - logs Liters/minute every minute
 - automatically enabled during sampling

Active Flow Controller

- Developed by APL with MBARI
 - combines digital flow meter with pump speed controller
 - ESP sets desired input loop flow rate during sampling
 - Powered by External Sampler power switch
 - Communicates on External Sampler serial port

Active Flow Controller Config

- Configuration in configure.rb:

```
class Sampler
  const :PumpPrimeDelay, 3.minutes #prime before sampling
  const :PrimeFlow,      2.0..7.0   #Liters/min flow range while priming
  const :MaintainFlow,  1.0..3.5   #flow range while sampling
end
```

- Configuration in corecal.rb:

```
flowScale = 0.00276964 #scales counts/min flow to liters/min
@@flow = Skew flowScale,-flowScale, format: '%.1fL/min', threshold:0,
        driver: Gauge::Flow::UWFlowControl,
        port: '/dev/extSampler', switch: :sampler
```


Active Flow Controller Use

- See lib/gauge/uwflowcontrol.rb

-> `Flow.start` #start monitoring flow

-> `Flow.selectIntake :shallow` #select shallow intake

-> `Flow.selectIntake :deep` #select deep intake

-> `Flow.intake` #verify intake selection

:deep

-> `Flow.prime` #run pump at faster priming flow rate

-> `Flow.maintain` #run pump at slower sampling flow rate

-> `Flow.rate` #returns current flow rate in Liters/minute

3.2

-> `Flow.disablePump` #stop pump

-> `Flow.stop` #stop pump and stop monitoring flow

